

# Mapping read data

MAQ, BOWTIE2, BWA, SOAP2, SSAHA2, ELAND2,  
RMAP, NovoAlign, Mosaik

Mapping read data :

Bowtie2

<http://bowtie-bio.sourceforge.net/bowtie2/index.shtml>

<http://bowtie-bio.sourceforge.net/bowtie2/manual.shtml>

# Building genome index using Bowtie2

► `bowtie2-build [options] <reference_in> <bt2_base>`

Command	What does it mean?
<code>&lt;reference_in&gt;</code>	A comma-separated list of files containing reference sequences
<code>&lt;bt2_base&gt;</code>	The basename of the index file to write (would be written in .bt2 format)

# Using Bowtie2

► bowtie2 [options] -x <bt2-idx> {-1 <m1> -2 <m2> | -U <r> | --interleaved <i>} -S [<sam>]

Command	What does it mean?
-x <bt2-idx>	Basename for Index file for the reference genome
-1 <m1>	File containing Paired end read 1 / mate 1 sequences
-2 <m2>	File containing Paired end read 2 / mate 2 sequences
-U <r>	File containing Unpaired read sequences
-- interleaved <i>	File containing read 1 and read 2 sequences interleaved
-S <sam>	Output file containing alignments in .sam format
-q	Reads in <m1>, <m2> , <r> are .fastq files

# Bowtie2 : Input options

Command	What does it mean?
-5/--trim5 <int>	Trim <int> bases from 5' (left) end of each read before alignment (default: 0).
-3/--trim3 <int>	Trim <int> bases from 3' (right) end of each read before alignment (default: 0).
--trim-to [3: 5:]<int>	Trim reads exceeding <int> bases. Bases will be trimmed from either the 3' (right) or 5' (left) end of the read. If the read end is not specified, bowtie 2 will default to trimming from the 3' (right) end of the read. <u>--trim-to</u> and <u>-3/-5</u> are mutually exclusive.
--phred33	Input qualities are ASCII chars equal to the <u>Phred quality</u> plus 33. This is also called the “Phred+33” encoding, which is used by the very latest Illumina pipelines.
--phred64	Input qualities are ASCII chars equal to the <u>Phred quality</u> plus 64. This is also called the “Phred+64” encoding.

# Bowtie2 : Run options

Command	What does it mean?
--end-to-end mode	Alignment from one end of the read to the other end
--very-fast	
--fast	
--sensitive	(default in --end-to-end mode)
--very-sensitive	
--local mode	Not entire read is required to align (soft clipping/trimming)
--very-fast-local	
--fast-local	
--sensitive-local	(default in --local mode)
--very-sensitive-local	

# Bowtie2 : Output options

Command	What does it mean?
<code>--fr/--rf/--ff</code>	The upstream/downstream mate orientations for a valid paired-end alignment against the forward reference strand. E.g., if <code>--fr</code> is specified and there is a candidate paired-end alignment where mate 1 appears upstream of the reverse complement of mate 2 Default: <code>--fr</code> (appropriate for Illumina's Paired-end Sequencing Assay).
<code>--no-mixed</code>	By default, when bowtie2 cannot find a concordant or discordant alignment for a pair, it then tries to find alignments for the individual mates. This option disables that behavior.
<code>--no-discordant</code>	By default, bowtie2 looks for discordant alignments if it cannot find any concordant alignments. A discordant alignment is an alignment where both mates align uniquely, but that does not satisfy the paired-end constraints ( <code>--fr/--rf/--ff</code> ). This option disables that behavior.

# Bowtie2 : Output options

Command	What does it mean?
<code>-l/--minins &lt;int&gt;</code>	The minimum fragment length for valid paired-end alignments. E.g. if <code>-l 60</code> is specified and a paired-end alignment consists of two 20-bp alignments in the appropriate orientation with a 20-bp gap between them, that alignment is considered valid (as long as <code>-X</code> is also satisfied).
<code>-X/--maxins &lt;int&gt;</code>	The maximum fragment length for valid paired-end alignments. E.g. if <code>-X 100</code> is specified and a paired-end alignment consists of two 20-bp alignments in the proper orientation with a 60-bp gap between them, that alignment is considered valid (as long as <code>-l</code> is also satisfied). Default: 500.