# Data Formats Data Quality Read Mapping

Riddhiman Dhar Department of Biotechnology IIT Kharagpur

#### Data formats

- Roche 454
  - FASTA format

- ► Illumina
  - FASTQ format

#### Data formats

- PacBio
  - basecall HDF5 format earlier
  - currently BAM format (convertible to fastq format)
- Ion-torrent
  - DAT files (Raw data)
  - BAM, FASTQ or VCF format
- Nanopore
  - HDF5 format
  - FAST5 format

### Fastq format

► Fasta format with Quality scores

#### Fastq format

- ► Fasta format with Quality scores
- @D00733:181:CAH6EANXX:8:2210:1499:2056 1:N:0:GTAGAG
- NCTTTGTACTATGACCGATACACTCAACCGGCGAAAAGTGGAACTTGAGAATTGATGTCTTCATCTTATT ATCTGTCTCTTATACACATCTCCGAGCCCACGAGACGTAGAGGGAATCTCGTATGC
- +
- @D00733:181:CAH6EANXX:8:2210:1668:2177 1:N:0:GTAGAG
- +

#### Header: @D00733:181:CAH6EANXX:8:2210:1668:2177 1:N:0:GTAGAG

Element	Requirements	Description			
@	@	Each sequence identifier line starts with @			
<instrument></instrument>	Characters allowed: a–z, A–Z, 0–9 and underscore	Instrument ID			
<run number=""></run>	Numerical	Run number on instrument			
<flowcell id=""></flowcell>	Characters allowed: a–z, A–Z, 0–9				
<lane></lane>	Numerical	Lane number			
<tile></tile>	Numerical	Tile number			
<x_pos></x_pos>	Numerical	X coordinate of cluster			
<y_pos></y_pos>	Numerical	Y coordinate of cluster			
<read></read>	Numerical	Read number. 1 can be single read or Read 2 of paired-end			
<is filtered=""></is>	Y or N	Y if the read is filtered (did not pass), N otherwise			
<control number=""></control>	Numerical	0 when none of the control bits are on, otherwise it is an even number.			
		On HiSeq X and NextSeq systems, control specification is not performed and this number is always 0.			
<sample number=""></sample>	Numerical	Sample number from sample sheet			

## Quality score - ASCII table

Dec	c Char		Dec	Char		Char		Char
0	NUL	(null)	32	SPACE	64		96	
1	SOH	(start of heading)	33	!	65	A	97	a
2		(start of text)	34	"	66	В	98	b
3	ETX	(end of text)	35	#	67	С	99	С
4	EOT	(end of transmission)		Ş	68	D	100	d
5	ENQ	(enquiry)	37	%	69	E	101	е
6	ACK	(acknowledge)	38	&	70	F	102	f
7	BEL	(bell)	39	•	71	G	103	g
8	BS	(backspace)	40	(	72	H	104	h
9	TAB	(horizontal tab)	41	)	73	I	105	i
10	LF	(NL line feed, new line)	42	*	74	J	106	j
11	VT	(vertical tab)	43	+	75	K	107	k
12	FF	(NP form feed, new page)	44	,	76	L	108	1
13	CR	(carriage return)	45	_	77	M	109	m
14	SO	(shift out)	46		78	N	110	n
15	SI	(shift in)	47	/	79	0	111	0
16	DLE	(data link escape)	48	0	80	P	112	р
17			49	1	81	Q	113	q
18	DC2	(device control 2)	50	2	82	R	114	r
19	DC3	(device control 3)	51	3	83	S	115	S
20	DC4	(device control 4)	52	4	84	T	116	t
21	NAK	(negative acknowledge)	53	5	85	U	117	u
22	SYN	(synchronous idle)	54	6	86	V	118	v
23		(end of trans. block)	55		87	W	119	W
24	CAN	(cancel)	56	8	88	X	120	X
25	EM	(end of medium)	57	9	89	Y	121	У
26	SUB	(substitute)	58	:	90	Z	122	Z
27	ESC	(escape)	59	;	91	[	123	{
28	FS	(file separator)	60	<	92	\	124	
29	GS	(group separator)	61	=	93	]	125	}
30	RS	(record separator)	62	>	94	^	126	~
31	US	(unit separator)	63	?	95	_	127	DEL

#### HDF5 format

- ► Hierarchical Data Format (HDF)
- ► Nested data
  - Groups
  - Datasets
  - Attributes

#### FAST5 format

- Specifically structured HDF5 data
- "Raw" data and "Analysis" data
- Raw data values of pico-amp currents in nanopore
- ► Analysis data basecall data


## Checking NGS Data Quality

#### Checking NGS Data Quality


#### FastQC

https://www.bioinformatics.babraham.ac.uk/projects/fastqc/


#### Per base quality score



#### Per base quality score




#### Per base quality score




Bad data!


## Per sequence quality scores




## Per sequence quality scores




## Per sequence quality scores




#### Per base sequence content




## Per sequence GC content




#### Per Base N content




#### Sequence length distribution




## **Duplicate Sequences**



#### Over-represented k-mers



#### Adapter content



#### Read trimming tools

Trimmomatic

http://www.usadellab.org/cms/?page=trimmomatic

bbduk

https://jgi.doe.gov/data-and-tools/bbtools/bb-tools-userguide/bbduk-guide/

Cutadapt

https://pypi.org/project/cutadapt/1.3/